


Infection

Haemolytic uraemic syndrome caused by Shiga toxin-producing *E. coli*



Clinical information and laboratory results

A two-year-old patient presented at the hospital with watery stool, which later became bloody with fever appearing.

The blood count showed an HGB of 6.8 g/dL (4.22 mmol/L), resulting in an 'Anaemia' flag on the analyser and a PLT-F count of $27 \times 10^9/L$ triggering the 'Thrombocytopenia' flag. On day 4 after admission, schistocytes were present at 4.0% on the blood smear.

Testing results further showed acute kidney injury and elevated lactate dehydrogenase (LD), as well as an increase in WBC and C-reactive protein (CRP). This led to the suspicion of HUS. Antibodies to *Escherichia coli* O157 lipopolysaccharide were found, confirming the diagnosis of haemolytic uraemic syndrome (HUS) due to Shiga toxin-producing *E. coli*.




Fig. 1 Red blood cells classified as schistocytes by the Advanced RBC Application of the DI-60, an automated digital imaging analyser. The blood sample was obtained and analysed on day 5 after admission.

Result interpretation

Escherichia coli strain O157 is a type of enterohaemorrhagic *E. coli* (EHEC) [1] that is known to cause haemolytic uraemic syndrome (HUS). It colonises the gastrointestinal tract and produces a cytotoxin called Shiga toxin or Stx [2]. Stx can lead to a remodelling of cell expression. This can result in activation of platelets [3], expression of adhesion molecules [4] and inflammatory chemokines [5]. The cytotoxic effect of Stx is potentiated, promoting the adhesion of WBC to endothelial cells, causing thrombosis and tissue damage. Thrombosis in the capillaries subsequently leads to mechanical haemolysis [6], so patients present with fragmented RBC as well as decreased platelet counts.




Fig. 2 The patient's PLT histogram showing interferences at the upper discriminator on day 5 after admission.




Fig. 3 RET scattergram with an RBC population (blue) extended towards the low FSC area, indicating fragmented RBC.

References

- [1] Gambushe SM et al. (2022): Review of *Escherichia coli* O157:H7 Prevalence, Pathogenicity, Heavy Metal and Antimicrobial Resistance, African Perspective. *Infect Drug Resist.* 15: 4645–4673.
- [2] Tarr PI et al. (1995): *Escherichia coli* O157:H7: Clinical, Diagnostic, and Epidemiological Aspects of Human Infection. *Clin Infect Dis.* 20(1): 1–10.
- [3] Karpman D et al. (2001): Platelet activation by Shiga toxin and circulatory factors as a pathogenetic mechanism in the hemolytic uremic syndrome. *Blood.* 97(10): 3100–3108.
- [4] Morigi M et al. (1995): Verotoxin-1 promotes leukocyte adhesion to cultured endothelial cells under physiologic flow conditions. *Blood.* 86(12): 4553–4558.
- [5] Zoja C et al. (2002): Shiga toxin-2 triggers endothelial leukocyte adhesion and transmigration via NF-kappaB dependent up-regulation of IL-8 and MCP-1. *Kidney Int.* 62(3): 846–856.
- [6] Bommer M et al. (2018): The Differential Diagnosis and Treatment of Thrombotic Microangiopathies. *Dtsch Arztebl Int.* 115(19): 327–334.