SEED Hematologie

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Literatura

- [1] **Haeckel R, Wosniok W and Arzideh F.** (2007): A plea for intra-laboratory reference limits. Part 1. General considerations and concepts for determination. *Clin Chem Lab Med* 45 : 1033–1042.
- [2] **Henny J** (2009): The IFCC recommendations for determining reference intervals: strengths and limitations. *LaboratoriumsMedizin* 33 : 45–51.
- [3] **Cerioti F, Hinzmann R and Panteghini M** (2009): Reference intervals: the way forward. *Ann Clin Biochem* 46 : 8–17.
- [4] **Wakeman L, Al-Ismail S, Benton A et al.** (2007): Robust, routine haematology reference ranges for healthy adults. *Int J Lab Hematol* 29 : 279–283.
- [5] **Heil W and Ehrhardt V** (2008): Reference ranges for adults and children. Pre-analytical considerations. Mannheim: Roche Diagnostics.
- [6] **Bain BJ and England JM** (1975): Normal haematological values: sex difference in neutrophil count. *Br Med J.* 5953 : 306–309.
- [7] **Solberg HE** (2004): The IFCC recommendation on estimation of reference intervals. The RefVal program. *Clin Chem Lab Med* 42 : 710–714.
- [8] **Hong J, Min Z, Bai-Shen P, Jie Z, Ming-Ting P, Xian-Zhang H, Xiao-Ke H, Lan-Lan W, Xin Z, Wei G, Rui Q, Wen-Xiang C, Xin-Zhong W, Yue-Yun M, Hong S** (2014): Investigation on Reference Intervals and Regional Differences of Platelet Indices in Healthy Chinese Han Adults. *J Clin Lab Anal*; early online.
- [9] **Qiao R, Yang S, Yao B, Wang H, Zhang J, Shang H** (2014): Complete blood count reference intervals and age- and sex-related trends of North China Han population. *Clin Chem Lab Med*; early online.
- [10] **Ambayya A, Su AT, Osman NH, Nik-Samsudin NR, Khalid K, Chang KM, Sathar J, Rajasuriar JS, Yegappan S** (2014): Haematological Reference Intervals in a Multiethnic Population. *PLoS One* 18; 9(3) : e91968.
- [11] **Sehgal KK, Tina D, Choksey U, Dalal RJ, Shanaz KJ** (2013): Reference range evaluation of complete blood count parameters with emphasis on newer research parameters on the complete blood count analyzer Sysmex XE-2100. *Indian J Pathol Microbiol* 56(2) : 120–4.
- [12] **El Graoui O, Ellaabi S, Igala M, Elaloussi FZ, Nourchafi N, Mifdal H, Serhier Z, Othmani MB, Mikou KA, Oukkache B** (2014): Hematology reference intervals in Moroccan population. *Clin Lab* 60(3) : 407–11.
- [13] **Pekelharing JM, Hauss O, de Jonge R, Lokhoff J, Sodikromo J, Spaans M, Brower R, de Lathouder S, Hinzmann R** (2010): Haematology reference intervals for established and novel parameters in healthy adults. *Sysmex Journal International* Vol. 20. No. 1.
- [14] **Sysmex Corporation, Japan** (2014): Reference ranges analysis document for XN series.
- [15] **Park SH, Park CJ, Lee BR, Kim MJ, Han MY, Cho YU, Jang S** (2015): Establishment of Age- and Gender-Specific Reference Ranges for 36 Routine and 57 Cell Population Data Items in a New Automated Blood Cell Analyzer, Sysmex XN-2000. *Ann Lab Med*.36(3) : 244–9.
- [16] **Imperiali CE, Arbiol-Roca A, Sanchez-Navarro L, Dastis-Arias M, Lopez-Delgado JC, Cortes-Bosch A, Sancho-Cerro A, Dot-Bach D** (2018): Reference interval for immature platelet fraction on Sysmex XN haematology analyser in adult population. *Biochem Med*.28(1) : 010708.
- [17] **Arbiol-Roca A, Imperiali CE, Montserrat MM, Cerro AS, Bosch de Basea AC, Navarro LS, Dot Bach D, Politi JV** (2018): Reference intervals for a complete blood count on an automated haematology analyser Sysmex XN in healthy adults from the southern metropolitan area of Barcelona. *EJIFCC*.29(1) : 48–54.
- [18] **Florin L, Oyaert M, Vandevenne M, Van den Bossche J, Vanden Driessche M, Claeys R, Marcelis L, Robbrecht J, Jacobs N, Nijs A, Janssen K, Van Ruymbeke A, de Mûelenaere G, Stove V** (2020): Establishment of common reference intervals for hematology parameters in adults, measured in a multicenter study on the Sysmex XN-series analyzer. *Int J Lab Hematol.* 42(3) : e110–e115.